

madalù

MADE IN ITALY

MADE IN ITALY



HIGH QUALITY



SHORT CHAIN



GUARANTEE AND
SAFETY PRODUCTS

WHEN NO ONE HAD
THOUGHT
WE HAD ALREADY
CREATED



10 YEARS OF KNOW HOW
IN AN INNOVATIVE
PRODUCT , THE FIRST TO
ENTER IN THE ITALIAN
MARKET

DIRECT PURCHASE FROM
MANUFACTURER

GREATER MARGIN AND
MORE COMPETITIVENESS
COMPARED TO SIMILAR
PRODUCTS



HOW WE PRODUCE

WHITE ROOM PRODUCTION

The production of the Mamicup® menstrual Cup takes place in a **white room**, 100% sterile, packaged one by one inside.
The guarantee of quality and certainty of safety of the product is in this way guaranteed 100% on every single piece.



Women want to be sure of the product they buy, especially in this area

Increased sales and product demand

Manufactured in FDA certified medical grade silicone of the highest quality

Designed by a woman, Mamicup® has been on the market for over 10 years and has been tested and approved by thousands of women all over the world



100%
MADE IN
ITALY

SAFETY
I WANT AN HIGH QUALITY
PRODUCT

- Hypoallergenic
- Made In Italy
- Latex Free
- **Ecological**
- Phthalates Free
- Practice
- Odorless
- **Reusable**
- Economic
- Comfortable



STOP to chemicals in contact with the vaginal mucosa
Mamicup® is manufactured in 100% medical grade platinum silicion.
It is also perfect for women who are allergic to latex or who have very sensitive skin.

STOP the dryness caused by the swab that absorbs the flow but also the protective mucus and alters the vaginal PH

STOP irritation, redness, discomfort, humidity



STOP to waste every month:
No more sanitary napkins, pads, applicators, sachets or plastic strips are thrown.

As many as 45,000 million external and internal sanitary napkins are thrown around the world every year.





STOP for every concern:
the Mamicup® is emptied only every 4/8 hours. Three times more autonomy than a traditional pad or pad
STOP for bulky packages:
Mamicup® is used day and night, throughout the cycle. When you travel, just pack a small bag


STOP to renunciations:
with Mamicup® you move in total freedom. Enter and exit the water without having to change. Play sports, you are active the way you like it. The cycle is no longer a limit.

STOP to emergencies:
wherever you are with Mamicup® you are calm even if the cycle suddenly arrives

COLORS AND SIZES

MAMICUP® IS PRODUCED IN 2 SIZES AND 4 COLORS

I have had children?	I have less than 30 years	I have more than 30 years
NO	"M"	"L"
YES with birth vaginally	"L"	"L"
YES by caesarean section	"M"	"L"




M (small)

capacity
25 ml

diameter: 40 mm

length of the stem:
23 mm




length (without stem): 50 mm

L (big)

capacity
35 ml

diameter: 45 mm

length of the stem:
18 mm




length (without stem): 57 mm

M

EVERY SIZES 4 COLORS

L



TRANSPARENT



LILAC




PINK



TURQUOISE

HOW WE USE

1



FOLD

Push down the top to open it


2



INSERT

Delicately put yourself in a comfortable position


3



EXTRACT

After 4, maximum 8 hours, gently pulling the stem

4



EMPTY

Rinse it and insert it again

How we use



COPPETTA MESTRUALE
ASSORBENTI LAVABILI
SPUGNE DI MARE NATURALI



SILVIA CARIOLATO

Business Development Manager

Via Piave 26/B - Fossalta di Portogruaro (VE)

silvia.cariolato@coppetta-mestruale.com

Uff. +39 0421 184 77 10

Cell. +39 389 5555 367

